remel

Pathotec® CO (Cytochrome Oxidase)

(English)

INTENDED USE

Remel Pathotec® CO (Cytochrome Oxidase) is a reagent impregnated paper strip recommended for use in qualitative procedures for the rapid detection of cytochrome oxidase activity in microorganisms.

SUMMARY AND EXPLANATION

The enzyme, cytochrome oxidase (CO), is produced in large amounts by some genera and species of bacteria. However, this enzyme is not produced by members of the family *Enterobacteriaceae*. Furthermore, typical genera and species of this family of microorganisms are not only cytochrome oxidase negative but also nitrate reductase positive. When Pathotec® CO is used with a nitrate reductase test, the results of these two tests generally indicate whether an isolate belongs to the family *Enterobacteriaceae* (exceptions are *Stenotrophomonas maltophilia*¹ and some members of the genus *Pasteurella*).

Pathotec® CO has been found to compare well with conventional procedures for detecting cytochrome oxidase activity in microorganisms and offers the following advantages: (1) Accuracy – complete agreement between Pathotec® CO and conventional methods^{2,3,4,5}; (2) Speed – following primary isolation, Pathotec® CO can be read 30 seconds after inoculation; (3) Flexibility – blends conveniently with current procedures and requires only a small storage area compared to conventional media; (4) Reliability – reagents which duplicate conventional procedures are applied to paper, dried, and packaged in desiccant-capped, amber vials to ensure stability.

PRINCIPLE

The enzyme cytochrome oxidase couples with α -naphthol and N,N-dimethyl-p-phenylenediamine in the presence of oxygen to form a bright blue product, indophenol blue, and water.

REAGENTS

Reactive ingredients: N,N-dimethyl-p-phenylenediamine
 α -naphthol

PRECAUTIONS

This product is For *In Vitro* Diagnostic Use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, and media after use. Directions should be read and followed closely.

STORAGE

This product is ready for use and no further preparation is necessary. The product should be stored in its original container at 2-8°C until used. Allow product to come to room temperature before use. Protect product from light.

PRODUCT DETERIORATION

The product should not be used if (1) the color has changed, (2) the expiration date has passed, or (3) there are other signs of deterioration. Protect strips from moisture, heat, and light. Promptly replace the cap and store at 2-8°C.

SPECIMEN COLLECTION, STORAGE, AND TRANSPORTATION

Specimens should be collected and handled following recommended procedures.^{6,7}

MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms, (6) Forceps, (7) Platinum loop or needle, or wooden applicator stick.

PROCEDURE

Clinical specimens should be inoculated onto primary isolation media as soon as possible after receipt. Specimens must be carefully streaked to obtain well-defined, isolated colonies.^{1,6}

NOTE: Colony sizes of *Enterobacteriaceae* on primary media vary substantially. One medium-sized colony will contain a sufficient number of cells for the Pathotec® CO test.

1. Place Pathotec® CO Test Strip on a paper towel or in a petri dish.
2. Select a well-isolated colony from a primary isolation medium.
3. Using a platinum loop or needle, or wooden applicator stick, smear the growth of an 18-24 hour old culture onto the pink-gray reagent band of the Pathotec® CO strip.
4. Allow to stand approximately 30 seconds.
5. Observe the color of the inoculated area.

INTERPRETATION OF THE TEST

Positive Test - Blue color development

Negative Test - Absence of blue color

QUALITY CONTROL

All lot numbers of Pathotec® CO have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

CONTROL	INCUBATION	RESULTS
<i>Pseudomonas aeruginosa</i> ATCC® 27853	Aerobic, 30 seconds @ Room Temperature	Positive
<i>Alcaligenes faecalis</i> ATCC® 8750	Aerobic, 30 seconds @ Room Temperature	Positive
<i>Escherichia coli</i> ATCC® 25922	Aerobic, 30 seconds @ Room Temperature	Negative

PERFORMANCE CHARACTERISTICS

The performance of the Pathotec® CO has been confirmed by several investigators. Complete agreement has been reported between Pathotec® CO and conventional procedures.^{2,3,4,5}

LIMITATIONS

1. Use a fresh culture (18-24 hours old). Older cultures or cultures grown on old or deteriorated media will have decreased levels of enzyme activity and may give false negative test results. Direct use of broth cultures is not recommended.⁸
2. Select colonies for testing which have not been cultured on selective media or media containing glucose, as selective agents and fermentation may inhibit oxidase enzyme activity resulting in possible false-negative results.⁸
3. A platinum inoculating loop or needle, or a wooden applicator stick is recommended for use in oxidase testing. The presence of any trace of iron (nichrome) can catalyze an oxidase reaction, resulting in a false-positive reaction.⁸

BIBLIOGRAPHY

1. Edwards, P.R. and W.H. Ewing. 1972. Identification of *Enterobacteriaceae*. 3rd ed. Burgess Publishing Co., Minneapolis, MN.
2. Blazevic, D.J., P.C. Schreckenberger, and J.M. Matsen. 1973. Appl. Microbiol. 26:886-889.
3. Edberg, S.C., M. Novak, H. Slater, and J.M. Singer. 1975. J. Clin. Microbiol. 2:469-473.
4. Rosner, R. 1973. Appl. Microbiol. 26:890-893.
5. Smith, P.B., D.L. Rhoden, and K.M. Tomfohrde. 1975. J. Clin. Microbiol. 1:359-362.
6. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Pfaller, and R.H. Tenen. 2003. Manual of Clinical Microbiology. 8th ed. ASM, Washington, D.C.
7. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2002. Bailey and Scott's Diagnostic Microbiology. 11th ed. Mosby, St. Louis, MO.
8. MacFaddin, J.F. 2000. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Lippincott Williams & Wilkins, Philadelphia, PA.

PACKAGING

REF 38191, Pathotec® CO 100 Strips/Pk (4 x 25/Vial)

Symbol Legend

REF	Catalog Number
IVD	In Vitro Diagnostic Medical Device
LAB	For Laboratory Use
	Consult Instructions for Use (IFU)
	Temperature Limitation (Storage Temp.)
LOT	Batch Code (Lot Number)
	Use By (Expiration Date)
EC REP	European Authorized Representative

Pathotec® is a registered trademark of Remel Inc.
ATCC® is a registered trademark of American Type Culture Collection.

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Pathotec® es una marca registrada de Remel Inc.